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Stimulation suit / Stimulation vest / Strap electrodes

StimaWELL® EMS

This instruction for use is valid for all stimulation suits, stimulation vests and strap electrodes of Pierenkemper GmbH.

Stimulation vests (12 channel):

462319, 462318, 462317, 462316, 462315, 462314, 462313, 462312, 462311, 462310

Stimulation vests (6 channel):

462329, 462328, 462327, 462326, 462325, 462324, 462323, 462322, 462321, 462320

Stimulation suits (12 channel):

462220, 462221, 462222, 462223, 462224, 462225, 462226, 462227, 462290, 462291, 462354, 462355, 462356, 462357, 462358

Stimulation suits (6 channel):

462191, 462190, 462189, 462188, 462187, 462186, 462185, 462184, 462183, 462178

Strap electrodes:

462292, 462293, 462294, 462252, 462253, 462254, 462193, 462194, 462195, 462267, 462268, 462269

General information

Please note, that under the term „stimulation suit“ we include the stimulation vest.

CONTENT

Use and safety Information	13
Recommendation for commercial use	15
Washing instructions, cleaning and care	16
Warranty/Repair.....	17
Material wearing.....	17
Step by step instructions	18

USE AND SAFETY INFORMATION

For your own safety or the safety of the user, please follow the instructions for use and safety information below. A trouble-free EMS training can only be guaranteed if you follow these instructions. Not following these instructions may impair the functionality of the products, the success of the training and may result in damage to the device or lead to health hazards. If the following instructions are not observed, we exclude quality claims on the product. Please also observe the instructions for use of your electrostimulation device. Use the stimulation suit/strap electrodes in the following order: (In addition, after reading the following points, you can find an illustrated summary at the bottom of page 18/19).

Getting started

- » The stimulation suit consists of polyamide and elastane, the conductive surfaces of the suit are made of nylon and silver. If you are allergic to any of the ingredients, use the product with caution.
- » Please wash suits/strap electrodes according to the washing instructions (see below) before first use.
- » For commercial use: Before each operation, especially after washing, check the conductivity of the electrodes. Insufficient conductivity hinders the success of the training and can lead to health hazards. A multimeter can be used for checking the conductivity of the electrodes. Place the two sensors of the multimeter at a maximum distance from each other on the electrode surface - if possible at a distance of 10 cm. If the resistance is higher than 30 Ohm, the electrode is not conductive enough and must be replaced. You can purchase a multimeter at www.stimawell-ems.de/en
- » Always follow the care and cleaning instructions, otherwise the quality and functionality of the products cannot be guaranteed.
- » First, put on the EMS underwear. Stimulation without underwear is also possible.
- » Before putting the suit on, fasten the distribution board onto the appropriate part of the suit for easy handling.

Moistening of the electrodes (silver-coated surfaces)

- » Moisten the electrodes before slipping over and tightening the stimulation suit or the strap electrodes. Moisten only the electrodes (not the entire suit but only the electrode surfaces). Moisten the electrodes in the correct ratio. Electrodes do not have to be dripping wet with the use of EMS underwear but may need to be moistened with a little more water as the water passes through the garment. For training without underwear, moistening with a wet sponge is sufficient. In case of hairiness you should moisturize a little more. Do not spray water neither on the EMS underwear nor on the suit fabric. Only spray the electrode surfaces in the suit.

- » Moisten the electrodes well and sufficiently with water before each workout to ensure optimum conductivity over the entire surface. Otherwise the functionality of the products cannot be guaranteed.
- » Do not spray the suit from the outside even if it is already worn. The whole suit gets wet and the current creeps through the whole suit fabric without causing proper muscle contractions.
- » Do not immerse the suit in water! Remember that electricity always takes the shortest route. If the whole suit is wet, it can happen that the current jumps from one electrode to the next and therefore no or too much stimulation takes place.
- » Particularly pay attention not to spray too much on the abdominal electrodes. Since the distance between these electrodes is small, the current could jump from one electrode to the other. This happens especially when the surrounding suit fabric or underwear is too wet.
- » Give the electrodes time to completely absorb the water. With completely new or freshly washed electrodes, the water can initially roll off the electrode material. Then slightly rub the water onto the electrodes with a sponge.
- » Only use original EMS underwear from StimaWELL. This is intentionally cut so short that the waistbands of the shorts and the shirt do not lie on top of each other. If using EMS underwear from other manufacturers or incorrectly selected sizes, these cuffs may overlap. The current does not find its way to the skin and no stimulation is perceived.
- » **Attention!** Inadequate or improper moistening may result in hazards.

Putting on the suits/strap electrodes

- » Carefully get in and out of the stimulation suit. If possible, get help getting in and out of the suit to avoid cable cracks! Do not step into the suit with shoes to avoid cable and seam cracks.
- » If the electrodes of the stimulation suit are not wet enough, open the zipper of the suit, if necessary partially pull out the suit - and then spray the electrodes on the inside by raising the suit.
- » The stimulation suit should always worn tightly to the body. However, it should be noted that excessive stretching can damage the material and cause a loss of elasticity of the suits.
- » After putting on the suit/strap electrodes, connect the strap electrodes to the suit using the connecting cables.
- » We recommend wearing an abdominal belt during EMS training to ensure better contact pressure of the abdominal and back electrodes. Wearing an abdominal belt is particularly important during dynamic workout, if you have a hollow back or with heavy body hair.
- » Now connect the distributor connecting cable to the distribution board.

After use

- » Before removing the suit, the first step is to disconnect the connection between the distribution board and the distributor connection cable.
- » Hang up the suit after each workout and let it dry completely. This is absolutely necessary to keep full future functionality of the electrodes. Complete drying ensures that no creep current can occur in the repeatedly moistened stimulation suit.
- » Don't wring out the suit. Internal cables could break as a result.
- » Defective suits should always be washed and returned to us with a description of the failure. A document for failure description is available on request.
- » Do not spray the electrode material with disinfectant spray or odour remover. This will destroy the electrode material!

Current irregularities

If you detect current irregularities during EMS training, follow the following instructions.

- » Measure the resistance on the push-buttons on the suit of the corresponding electrode (slots for distribution board on the suit). The electrode assignment plan is available on request. If the resistance of the push-button is higher than 5 Ohm, one or several internal wires are torn. The suit must be washed and returned to us with a description of the failure.
- » Check the push-button of the connecting cables (connection between the suit and the strap electrodes).
- » Check the correct fit of the stimulation suit and the strap electrodes. The electrodes should be positioned directly above the corresponding muscles. Abdominal or buttock electrodes which are too high or too low indicate an incorrect suit size.

RECOMMENDATION FOR COMMERCIAL USE

For a trouble-free commercial use you need 15 suits in different sizes, a complete set of strap electrodes (6 x narrow, 6 x wide strap electrode sets) as well as spare connecting cable for the connection suit/strap electrode and a spare distribution board. A studio operation with less than 15 suits and spare parts can only be carried out under the acceptance of possible disturbances and/or considerable delays, as each suit must be washed and left to dry after each use. We ask for your understanding that we cannot assume any liability in the case of non-observance.

WASHING INSTRUCTIONS, CLEANING AND CARE

Shades, colour differences and batik-like discolourations are not a manufacturing defect and do not represent a quality deficiency. To ensure that the conductivity is maintained even after several washing cycles, please always observe the following care instructions:

Before washing it is essential to close the Velcro fasteners on the arm and leg electrodes.

Washing cycle at maximum temperature 30°C. Only use StimaWELL EMS or Merusal detergent. StimaWELL EMS detergent can be purchased online at stimawell-ems.com. Do not wash at higher temperatures.

With the use of a normal household washing machine at 30°C more than 100 washing cycles are possible without any meaningful loss of electrode conductivity. 4 caps of Merusal EMS detergent should be used for a full washing machine load of 8 to 10 electrode suits with a pre-wash cycle at temperature 20° to 30°C. If the washing machine provides a soaking function then soak in the StimaWELL EMS or Merusal solution for 2 hours. If no soaking program is available then please use 8 caps of Merusal EMS detergent for a full machine load at a maximum temperature of 30°C and wash for at least 2 hours. Please do not use any other washing detergents, disinfection agents, chlorine based cleaning agents, softeners or textile fresheners.

No bleach. All bleaching agents will cause damage to the product. In no circumstances should chlorine be used because this will damage the silver fabric cladding and lead to a loss of function.

Do not iron.

Chemical cleaning is not allowed.

Dry after washing in the dryer with the program 'cupboard dry'.

Disinfection:

Do not use additional disinfection!. All types of disinfection will cause damage to the product – especially where chemical agents are involved. The silver fabric cladding acts anti-bacterially and therefore does not require prior additional disinfection.

Please note that a hygiene detergent (StimaWELL or Merusal EMS detergent) is not the same as a disinfectant detergent! Due to the unique surface properties we decline quality claims of incorrectly washed electrodes.

WARRANTY/REPAIR

The statutory warranty period of two years applies. The manufacturer and the dealer are not liable for damage caused by improper handling, non-compliance with the instructions for use or misuse.

The electrodes include the stimulation suits and the band electrodes (arm and leg electrodes).

The following applies to damage that is not covered by the warranty:

- » Electrodes can be repaired for a charge. The respective costs will be determined and communicated before the repair is carried out.
- » Electrodes sent for repair must be washed, regardless of whether it is a warranty case or not. Unwashed electrodes are sent back unprocessed.

MATERIAL WEARING

The electrodes and push buttons are consumables subject to natural wear. The durability, however, depends very much on the handling and especially on the correct care. Please follow the care instructions and only use StimaWELL or Merusal special detergent as a detergent. For these reasons, no warranty of durability can be given.

STEP BY STEP INSTRUCTIONS

1 Stimulation suit

Attach the distributor.

2 Stimulation suit

Check the connections.

3 Stimulation suit


Leg connections right and left: long cable on black push button. Short cable on silver push button.

4 Stimulation suit

Arm connections right and left: long cable on silver push button. Short cable on black push button.

5 Stimulation suit

Turn the suit inside out.

6 Stimulation suit

Moisten all electrodes with a sponge (also the side electrodes).

7 Stimulation suit

Turn the suit back to normal side, pull out the cable.

8 Stimulation suit

Carefully put on the suit without shoes (be careful with long nails).

9 Stimulation suit

Put on the first arm completely, then put on the other arm.

STEP BY STEP INSTRUCTIONS

10 Stimulation suit

Check that the electrodes are seated without wrinkles.

11 Stimulation suit

Hold zipper together and close (if necessary place electrodes correctly and check again the connections of the distributor).

12 Strap electrodes

Adjust each strap electrodes to the correct length, turn inside out and moisten.

13 Arm electrodes

Slip over upper arm electrodes and connect to short cable.

14 Arm electrodes

Slip over forearm electrodes and connect to long cable.

15 Leg electrodes

Slip over thigh electrodes and connect to short cable.

16 Leg electrodes

Slip over calf electrodes and connect to long cable.

17 Abdominal belt

Put on abdominal belt and adjust the distributor position.

18

Connect distributor and connection cable. Turn only slightly until the fastener engages (Caution! Pay attention to the colours: blue to blue and grey to grey).

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